

Anti-Drap1 Polyclonal Antibody For Wb And Elisa - Q14919

Item Number: CM0024160



Introduction

DRAP1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat DRAP1, a transcriptional corepressor. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	DRAP1 Rabbit pAb
Remarks/Alias	NC2-alpha; DRAP1
Species	Human
GeneID (Human)	10589
GeneID	10589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human DRAP1 (NP_006433.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14919
Protein Weight	22kDa
Shipping	Ice bag