

Anti-Drd1 Rabbit Polyclonal Antibody For Wb And Elisa - P21728

Item Number: CM0028119



Introduction

Rabbit polyclonal antibody against human DRD1 (Dopamine D1 receptor), suitable for Western blot and ELISA, with cross-reactivity to rat. Ideal for neuroscience and GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	DRD1 Rabbit pAb
Remarks/Alias	DADR; DRD1A; DRD1
Species	Human
GeneID (Human)	1812
GeneID	1812
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 338-446 of human DRD1 (NP_000785.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P21728
Protein Weight	49kDa
Shipping	Ice bag