


## Anti-Drd4 Polyclonal Antibody For Wb, Ihc-P, Elisa - P21917

Item Number: CM0013637



### Introduction

Anti-DRD4 polyclonal antibody raised in rabbit against a synthetic peptide within amino acids 100-200 of human DRD4. Validated for WB, IHC-P, ELISA, and cross-reacts with mouse and rat. Ideal for dopamine D4 receptor studies in neuroscience and behavior research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | DRD4 Rabbit pAb                                                                                                           |
| Remarks/Alias    | D4DR; DRD4                                                                                                                |
| Species          | Human                                                                                                                     |
| Gene ID (Human)  | 1815                                                                                                                      |
| Gene ID          | 1815                                                                                                                      |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human DRD4 (NP_000788.2). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IHC-P, ELISA                                                                                                          |
| Cross Reactivity | Mouse, Rat                                                                                                                |
| SWISS            | P21917                                                                                                                    |
| Protein Weight   | 44kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |