

Anti-Dreb1B Rabbit Polyclonal Antibody For Wb And Elisa - P93835

Item Number: CM0020388

Introduction

Rabbit polyclonal antibody targeting Arabidopsis thaliana DREB1B (CBF1), suitable for WB and ELISA. Ideal for cold stress transcription factor and freezing tolerance studies.

[Learn More](#)

Parameter	Specification
Product Name	DREB1B Rabbit pAb
Remarks/Alias	ATCBF1; C-repeat/DRE binding factor 1; DRE BINDING PROTEIN 1B; DREB1B; T30C3.11; TRANSCRIPTIONAL ACTIVATOR CBF1
Species	Arabidopsis thaliana
GeneID	828653
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of Arabidopsis thaliana DREB1B (NP_567721.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	P93835
Protein Weight	24kDa
Shipping	Ice bag