


Anti-Drg1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y295

Item Number: CM0026093



Introduction

Rabbit polyclonal antibody against human DRG1 (NEDD3), validated for WB and ELISA applications. Cross-reacts with mouse and rat orthologs. Ideal for research on GTP hydrolysis, microtubule polymerization, and mitotic spindle assembly.

[Learn More](#)

Parameter	Value
Product Name	DRG1 Rabbit pAb
Remarks/Alias	NEDD3; DRG1
Species	Human
Gene ID (Human)	4733
Gene ID	4733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 158-367 of human DRG1 (NP_004138.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y295
Protein Weight	41kDa
Shipping	Ice bag