

Anti-Drosha Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Nrr4

Item Number: CM0003561



Introduction

High-quality anti-DROSHA rabbit monoclonal antibody validated for Western blot, immunofluorescence, and ELISA. Ideal for miRNA biogenesis research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	DROSHA Rabbit mAb
Remarks/Alias	RN3; ETOH12; RNASEN; RANSE3L; RNASE3L; HSA242976; DROSHA
Species	Human
GeneID (Human)	29102
GeneID	29102
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DROSHA (Q9NRR4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NRR4
Protein Weight	159kDa
Shipping	Ice bag