


Anti-Drp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00429

Item Number: CM0027644



Introduction

Anti-DRP1 rabbit polyclonal antibody validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Detects human DRP1 and cross-reacts with mouse, rat. Ideal for studying mitochondrial fission, apoptosis, and peroxisomal division.

[Learn More](#)

Parameter	Value
Product Name	DRP1 Rabbit pAb
Remarks/Alias	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV
Species	Human
GeneID (Human)	10059
GeneID	10059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 411-710 of human DRP1 (NP_036193.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00429
Protein Weight	82kDa
Shipping	Ice bag