


Anti-Dsc2 Rabbit Polyclonal Antibody For Wb, Elisa - Q02487

Item Number: CM0016721



Introduction

High-quality DSC2 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Recognizes human desmocollin-2 (UniProt Q02487), a desmosome cadherin critical for cell adhesion and tissue integrity. Suitable for cardiac and intestinal research.

[Learn More](#)

Parameter	Value
Product Name	DSC2 Rabbit pAb
Remarks/Alias	DG2; DSC3; CDHF2; ARVD11; DGII/III; DSC2
Species	Human
Gene ID (Human)	1824
Gene ID	1824
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 520-690 of human DSC2 (NP_077740.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q02487
Protein Weight	100kDa
Shipping	Ice bag