

Anti-Dsg1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q02413

Item Number: CM0028018

Introduction

DSG1 Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence, and ELISA. Cross-reactive with human and rat DSG1 (~114 kDa). Ideal for desmosome adhesion, keratinocyte differentiation, and skin barrier studies.

[Learn More](#)

Parameter	Value
Product Name	DSG1 Rabbit pAb
Remarks/Alias	DG1; DSG; CDHF4; EPKHE; PPKS1; SPPK1; EPKHIA; DSG1
Species	Human
Gene ID (Human)	1828
Gene ID	1828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 800-1049 of human DSG1 (NP_001933.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q02413
Protein Weight	114kDa
Shipping	Ice bag