


Anti-Dsg2 Rabbit Polyclonal Antibody For Wb, Elisa - Q14126

Item Number: CM0023672



Introduction

Rabbit polyclonal antibody against human Desmoglein-2 (DSG2) validated for Western blot and ELISA. Ideal for cell adhesion, cardiac disease, and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	DSG2 Rabbit pAb
Remarks/Alias	HDGC, CDHF5, DSG2
Species	Human
GeneID (Human)	1829
GeneID	1829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 635-935 of human DSG2 (NP_001934.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14126
Protein Weight	122kDa
Shipping	Ice bag