

Anti-Dsg3/Desmoglein-3 Rabbit Polyclonal Antibody For IHC-P, Elisa - P32926

Item Number: CM0024535



Introduction

High-quality rabbit polyclonal antibody targeting human Desmoglein-3 (DSG3) for IHC-P and ELISA applications. Ideal for investigating desmosome adhesion, pemphigus vulgaris, and hair follicle biology.

[Learn More](#)

Parameter	Value
Product Name	DSG3/Desmoglein-3 Rabbit pAb
Remarks/Alias	PVA; ABOLM; CDHF6; Desmoglein 3
Species	Human
GeneID (Human)	1830
GeneID	1830
Immunogen	Recombinant protein corresponding to a sequence within amino acids 50-615 of human DSG3/Desmoglein-3 (NP_001935.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P32926
Protein Weight	108kDa
Shipping	Ice bag