

Anti-Dspp Polyclonal Antibody For Wb, Elisa - Q9Nzw4

Item Number: CM0025802

Introduction

Rabbit polyclonal antibody against human DSPP (Dentin sialophosphoprotein). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for dentin mineralization and biomineralization research.

[Learn More](#)

Parameter	Value
Product Name	DSPP Rabbit pAb
Remarks/Alias	DPP, DSP, DGI1, DMP3, DFNA39, DSPP
Species	Human
GeneID (Human)	1834
GeneID	1834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-146 of human DSPP (NP_055023.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZW4
Protein Weight	131kDa
Shipping	Ice bag