

Anti-Dspp Polyclonal Antibody For Wb - Q9Nzw4

Item Number: CM0023650

Introduction

Rabbit polyclonal antibody against human Dentin sialophosphoprotein (DSPP) for WB, IF-P, ELISA. Involved in dentinogenesis and biomineralization. Cross-reacts with human. Predicted molecular weight 131kDa.

[Learn More](#)

Parameter	Value
Product Name	DSPP Rabbit pAb
Remarks/Alias	DPP; DSP; DGI1; DMP3; DFNA39; DSPP
Species	Human
Gene ID	1834
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 16-230 of human DSPP (NP_055023.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	Q9NZW4
Protein Weight	131kDa
Shipping	Ice bag