

Anti-Dth8 Rabbit Polyclonal Antibody For Wb And Elisa - Q0J7P4

Item Number: CM0022864



Introduction

Rabbit polyclonal antibody against *Oryza sativa* DTH8 (NF-YB11), validated for WB and ELISA applications in rice; suitable for flowering time and growth regulation studies.

[Learn More](#)

Parameter	Value
Product Name	DTH8 Rabbit pAb
Remarks/Alias	Hd5; DTH8; GHD8; HAP3H; NFYB11; OsHAP3H; OsNF-YB11
Species	<i>Oryza sativa</i>
Gene ID	4344784
Immunogen	Recombinant protein of <i>Oryza sativa</i> DTH8
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q0J7P4
Protein Weight	30kDa
Shipping	Ice bag