

Anti-Dtna Polyclonal Antibody For Wb, Elisa - Q9Y4J8

Item Number: CM0014573



Introduction

Rabbit polyclonal antibody targeting human dystrobrevin alpha (DTNA). Validated for WB and ELISA, cross-reacts with human and mouse. Suitable for synapse and muscle biology research.

[Learn More](#)

Parameter	Value
Product Name	DTNA Rabbit pAb
Remarks/Alias	DTN; DRP3; DTN-A; LVNC1; D18S892E; DTNA
Species	Human
Gene ID (Human)	1837
Gene ID	1837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-480 of human DTNA (NP_001381.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y4J8
Protein Weight	84kDa
Shipping	Ice bag