

Anti-Dtnpb1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ev8

Item Number: CM0020893

Introduction

DTNBP1 rabbit polyclonal antibody validated for WB and ELISA in human and mouse. Dysbindin is a BLOC-1 complex component involved in lysosome-related organelle biogenesis and synaptic vesicle trafficking.

[Learn More](#)

Parameter	Value
Product Name	DTNBP1 Rabbit pAb
Remarks/Alias	SDY; DBND; HPS7; My031; BLOC1S8; DTNBP1
Species	Human
Gene ID	84062 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human DTNBP1 (NP_115498.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96EV8
Protein Weight	39kDa
Shipping	Ice bag