


Anti-Dtx1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Y01

Item Number: CM0021241



Introduction

Anti-DTX1 Rabbit Polyclonal Antibody (67kDa) for WB, ELISA. Validated in Human, Mouse. DTX1 is an E3 ubiquitin-protein ligase regulating Notch signaling and immune processes. Q86Y01.

[Learn More](#)

Parameter	Value
Product Name	DTX1 Rabbit pAb
Remarks/Alias	hDx-1; RNF140; DTX1
Species	Human
Gene ID (Human)	1840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human DTX1 (NP_004407.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86Y01
Protein Weight	67kDa
Shipping	Ice bag