

Anti-Dtx3L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Tdb6

Item Number: CM0020765



Introduction

Rabbit polyclonal antibody targeting human E3 ubiquitin-protein ligase DTX3L. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for DNA damage repair and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	DTX3L Rabbit pAb
Remarks/Alias	BBAP; RNF143; DTX3L
Species	Human
Gene ID (Human)	151636
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 601-740 of human DTX3L (NP_612144.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDB6
Protein Weight	84kDa
Shipping	Ice bag