


Anti-Dtymk Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P23919

Item Number: CM0031015



Introduction

Anti-DTYMK Rabbit Polyclonal Antibody validated for Western Blot, IF/ICC, and ELISA. Detects human, mouse, and rat thymidylate kinase, a key enzyme in dTTP biosynthesis. Suitable for nucleotide metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	DTYMK Rabbit pAb
Remarks/Alias	CDC8; TMPK; TYMK; CONPM; PP3731; DTYMK
Species	Human
GeneID (Human)	1841
GeneID	1841
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23919
Protein Weight	24kDa
Shipping	Ice bag