


Anti-Dusp13 Rabbit Polyclonal Antibody For Wb And Elisa - Q6B8I1/Q9Uii6

Item Number: CM0030948



Introduction

DUSP13 rabbit polyclonal antibody for WB and ELISA, targeting human, mouse, and rat DUSP13A/B isoforms. Supports research into phosphatase activity, ASK1-mediated apoptosis, and muscle-specific gene expression. Bulk supply available for IVD immunoassay development.

[Learn More](#)

Parameter	Value
Product Name	DUSP13 Rabbit pAb
Remarks/Alias	MDSP; TMDP; SKRP4; DUSP13
Species	Human
GeneID (Human)	51207
GeneID	51207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DUSP13 (NP_057448.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6B8I1/Q9Uii6
Protein Weight	6kDa/7kDa/9kDa/17kDa/20kDa/22kDa/27kDa
Shipping	Ice bag