


Anti-Dusp14 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95147

Item Number: CM0014018



Introduction

High-quality anti-DUSP14 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targeting dual specificity protein phosphatase 14 (MKP-6). Suitable for studying TCR signaling and MAP kinase inactivation.

[Learn More](#)

Parameter	Value
Product Name	DUSP14 Rabbit pAb
Remarks/Alias	MKP6; MKP-L; DUSP14
Species	Human
GeneID (Human)	11072
GeneID	11072
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DUSP14 (NP_008957.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95147
Protein Weight	22kDa
Shipping	Ice bag