

Anti-Dusp15 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H1R2

Item Number: CM0024412



Introduction

DUSP15 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA, reacts with human, mouse, rat. Explore dual-specificity phosphatase roles in signaling and myelination.

[Learn More](#)

Parameter	Value
Product Name	DUSP15 Rabbit pAb
Remarks/Alias	VHY; C20orf57; DUSP15
Species	Human
GeneID (Human)	128853
GeneID	128853
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-90 of human DUSP15 (NP_542178.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1R2
Protein Weight	32kDa
Shipping	Ice bag