

Anti-Dusp16 Rabbit Polyclonal Antibody For Wb, Elisa - Q9By84

Item Number: CM0013184



Introduction

Rabbit polyclonal antibody against DUSP16 (MKP-7), suitable for WB and ELISA, with reactivity in human, mouse, and rat. DUSP16 is a dual specificity phosphatase that inactivates MAP kinases. Ideal for research.

[Learn More](#)

Parameter	Specification
Product Name	DUSP16 Rabbit pAb
Remarks/Alias	MKP7, MKP-7, DUSP16
Species	Human
GeneID (human)	80824
GeneID	80824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 416-665 of human DUSP16 (NP_085143.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY84
Protein Weight	73kDa
Shipping	Ice bag