


Anti-Dusp22 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nrw4

Item Number: CM0025588



Introduction

Rabbit polyclonal antibody targeting human DUSP22 (JKAP/JSP-1), suitable for WB, IF/ICC, and ELISA. Explore JNK signaling and T-cell regulation with our high-quality antibody. Research use only.

[Learn More](#)

Parameter	Value
Product Name	DUSP22 Rabbit pAb
Remarks/Alias	VHX; JKAP; JSP1; MKPX; JSP-1; MKP-x; LMWDSP2; LMW-DSP2; DUSP22
Species	Human
GeneID (Human)	56940
GeneID	56940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human DUSP22 (NP_064570.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NRW4
Protein Weight	21kDa
Shipping	Ice bag