


Anti-Dusp5 Rabbit Polyclonal Antibody For Wb And Elisa - Q16690

Item Number: CM0016799



Introduction

Rabbit polyclonal antibody targeting human dual specificity protein phosphatase DUSP5 (hVH3). Suitable for Western blot and ELISA applications. Research use for ERK1 signaling and nuclear phosphatase studies. Molecular weight ~42 kDa.

[Learn More](#)

Parameter	Value
Product Name	DUSP5 Rabbit pAb
Remarks/Alias	DUSP; HVH3; DUSP5
Species	Human
GeneID (Human)	1847
GeneID	1847
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 315-384 of human DUSP5 (NP_004410.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16690
Protein Weight	42kDa
Shipping	Ice bag