


Anti-Dusp6 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16828

Item Number: CM0029903



Introduction

DUSP6 Rabbit pAb (Q16828) is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat DUSP6. Ideal for MAP kinase signaling and pain research studies.

[Learn More](#)

Parameter	Value
Product Name	DUSP6 Rabbit pAb
Remarks/Alias	HH19; MKP3; PYST1; DUSP6
Species	Human
GeneID (Human)	1848
GeneID	1848
Immunogen	Recombinant fusion protein containing a sequence within amino acids 302-381 of human DUSP6 (NP_001937.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16828
Protein Weight	42kDa
Shipping	Ice bag