


Anti-Dusp6 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q16828

Item Number: CM0003760



Introduction

High-quality Anti-DUSP6 Rabbit Monoclonal Antibody validated for WB, IP, IHC-P, and ELISA in human, mouse, rat. Targets DUSP6 (42 kDa), an ERK-specific MAP kinase phosphatase involved in cell differentiation and pain resolution.

[Learn More](#)

Parameter	Value
Product Name	DUSP6 Rabbit mAb
Remarks/Alias	HH19; MKP3; PYST1; DUSP6
Species	Human
GeneID (human)	1848
GeneID	1848
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 282-381 of human DUSP6 (Q16828).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16828
Protein Weight	42kDa
Shipping	Ice bag