


Anti-Dusp8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13202

Item Number: CM0021366



Introduction

High-quality rabbit polyclonal anti-DUSP8 antibody for WB, IF/ICC, and ELISA. DUSP8 is a dual specificity phosphatase that dephosphorylates MAP kinases, essential for signal transduction. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Specification
Product Name	DUSP8 Rabbit pAb
Remarks/Alias	HB5; HVH8; HVH-5; C11orf81; DUSP8
Species	Human
GeneID (human)	1850
GeneID	1850
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13202
Protein Weight	66kDa
Shipping	Ice bag