


Anti-Dusp9 (Mkp-4) Polyclonal Antibody For Wb, If/Icc, Elisa - Q99956

Item Number: CM0025603



Introduction

Rabbit polyclonal antibody targeting DUSP9 (MKP-4), validated for WB, IF/ICC, and ELISA with human and mouse reactivity. Ideal for ERK signaling research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	DUSP9 Rabbit pAb
Remarks/Alias	MKP4, MKP-4, DUSP9
Species	Human
Gene ID	1852
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DUSP9 (NP_001386.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99956
Protein Weight	42kDa
Shipping	Ice bag