

Anti-Dut Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P33316

Item Number: CM0029200



Introduction

Rabbit polyclonal DUT antibody for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Detects mitochondrial dUTPase to study nucleotide metabolism and DNA replication fidelity.

[Learn More](#)

Parameter	Value
Product Name	DUT Rabbit pAb
Remarks/Alias	BMFDMS; dUTPase; DUT
Species	Human
Gene ID (Human)	1854
Gene ID	1854
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 89-252 of human DUT (NP_001020419.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33316
Protein Weight	27kDa
Shipping	Ice bag