

Anti-Dvl1 Polyclonal Antibody For Wb, Elisa - O14640

Item Number: CM0023373



Introduction

Rabbit polyclonal antibody against human DVL1, valid for WB and ELISA, cross-reacts with mouse. Target involved in Wnt signaling; suitable for research use in developmental biology.

[Learn More](#)

Parameter	Value
Product Name	DVL1 Rabbit pAb
Remarks/Alias	DVL; DRS2; DVL1L1; DVL1P1; DVL1
Species	Human
GeneID (human)	1855
GeneID	1855
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 11-110 of human DVL1 (NP_001317240.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O14640
Protein Weight	75kDa
Shipping	Ice bag