

Anti-Dvl1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O14640

Item Number: CM0016444



Introduction

DVL1 Rabbit pAb is a polyclonal antibody targeting human DVL1 protein (UniProt O14640), suitable for WB, IHC-P, IF/ICC, and ELISA assays with cross-reactivity to mouse and rat. Ideal for Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	DVL1 Rabbit pAb
Remarks/Alias	DVL; DRS2; DVL1L1; DVL1P1; DVL1
Species	Human
GeneID (human)	1855
GeneID	1855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-640 of human DVL1 (NP_004412.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14640
Protein Weight	75kDa
Shipping	Ice bag