


Anti-Dvl2 Rabbit Polyclonal Antibody For If/Icc And Elisa - O14641

Item Number: CM0017524



Introduction

Rabbit polyclonal antibody against human DVL2, validated for IF/ICC and ELISA applications. Recombinant immunogen, cross-reacts with human, 79kDa protein. Ideal for studying canonical and non-canonical Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	DVL2 Rabbit pAb
Remarks/Alias	DVL2
Species	Human
GeneID (Human)	1856
GeneID	1856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human DVL2 (NP_004413.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O14641
Protein Weight	79kDa
Shipping	Ice bag