


Anti-Dvl3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92997

Item Number: CM0026355



Introduction

CamelBio's DVL3 rabbit polyclonal antibody targets human segment polarity protein dishevelled-3. Suitable for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for Wnt signaling studies.

[Learn More](#)

Parameter	Value
Product Name	DVL3 Rabbit pAb
Remarks/Alias	DRS3; DVL3
Species	Human
GeneID	1857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-650 of human DVL3 (NP_004414.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92997
Protein Weight	78kDa
Shipping	Ice bag