

Anti-Dynamin 3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uq16

Item Number: CM0015652



Introduction

Rabbit polyclonal antibody against human Dynamin 3 (Dyna III), validated for WB, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Ideal for vesicular trafficking and endocytosis research.

[Learn More](#)

Parameter	Value
Product Name	Dynamin 3 Rabbit pAb
Remarks/Alias	Dyna III; Dynamin 3
Species	Human
GeneID (human)	26052
GeneID	26052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-555 of human Dynamin 3 (NP_001265181.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UQ16
Protein Weight	98kDa
Shipping	Ice bag