

Anti-Dynamin 3 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Uq16

Item Number: CM0003369



Introduction

High-quality rabbit monoclonal antibody targeting human Dynamin-3 (DNM3). Validated for Western blot and ELISA applications. Detects human and mouse samples with predicted band at 98 kDa. Ideal for endocytosis and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	Dynamin 3 Rabbit mAb
Remarks/Alias	Dyna III; Dynamin 3
Species	Human
Gene ID (Human)	26052
Gene ID	26052
Immunogen	A synthetic peptide corresponding to amino acids 700-800 of human Dynamin 3 (Q9UQ16)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UQ16
Protein Weight	98kDa
Shipping	Ice bag