


Anti-Dync1H1 Rabbit Polyclonal Antibody For Wb, Elisa - Q14204

Item Number: CM0021628



Introduction

DYNC1H1 Rabbit Polyclonal Antibody validated for WB and ELISA. Specific to human DYNC1H1 (Cytoplasmic dynein 1 heavy chain 1, 532kDa), a motor protein for intracellular retrograde transport and mitotic spindle assembly. Ideal for cell biology and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	DYNC1H1 Rabbit pAb
Remarks/Alias	p22; DHC1; DNCL; DYHC; HL-3; CMT2O; DHC1a; DNCH1; DNECL; Dnchc1; CDCBM13; SMALED1; DYNC1H1
Species	Human
Gene ID (Human)	1778
Gene ID	1778
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human DYNC1H1 (NP_001367.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14204
Protein Weight	532kDa
Shipping	Ice bag