

Anti-Dync111 Polyclonal Antibody For Wb, Ihc-P, Elisa - O14576

Item Number: CM0024807

Introduction

Rabbit polyclonal antibody against human DYNC111, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for cytoplasmic dynein motor and intracellular transport studies. - O14576

[Learn More](#)

Parameter	Value
Product Name	DYNC111 Rabbit pAb
Remarks/Alias	DNC11; DNCIC1; DYNC111
Species	Human
GeneID (Human)	1780
GeneID	1780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human DYNC111 (NP_001129028.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14576
Protein Weight	73kDa
Shipping	Ice bag