


Anti-Dync112 Polyclonal Antibody For Wb, Elisa - Q13409

Item Number: CM0016846



Introduction

DYNC112 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA, reactive against human, mouse, and rat cytoplasmic dynein 1 intermediate chain 2, a key motor protein for intracellular transport.

[Learn More](#)

Parameter	Value
Product Name	DYNC112 Rabbit pAb
Remarks/Alias	IC2; DIC74; DNC12; NEDMIBA; DYNC112
Species	Human
Gene ID (Human)	1781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 352-612 of human DYNC112 (NP_001258717.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13409
Protein Weight	71kDa
Shipping	Ice bag