


Anti-Dync2H1 Polyclonal Antibody For Wb And Elisa - Q8Ncm8

Item Number: CM0014769



Introduction

High-quality rabbit anti-DYNC2H1 polyclonal antibody validated for WB and ELISA, reactive with human and rat. Sourced from recombinant immunogen (aa1-250). Ideal for cilia biogenesis and intraflagellar transport research.

[Learn More](#)

Parameter	Value
Product Name	DYNC2H1 Rabbit pAb
Remarks/Alias	ATD3; DHC2; DHC1b; DNCH2; DYH1B; SRTD3; SRP52B; hdhc11; DYNC2H1
Species	Human
Gene ID (Human)	79659
Gene ID	79659
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DYNC2H1 (NP_001073932.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8NCM8
Protein Weight	493kDa
Shipping	Ice bag