


Anti-Dynll1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P63167

Item Number: CM0013588



Introduction

High-quality rabbit polyclonal antibody targeting human DYNLL1 (dynein light chain 1). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for cytoskeleton and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	DYNLL1 Rabbit pAb
Remarks/Alias	LC8, PIN, DLC1, DLC8, LC8a, DNCL1, hdlc1, DNCLC1, DYNLL1
Species	Human
GeneID	8655
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-89 of human DYNLL1 (NP_003737.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63167
Protein Weight	10kDa
Shipping	Ice bag