


Anti-Dynll1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P63167

Item Number: CM0012808



Introduction

High-quality rabbit monoclonal antibody against human DYNLL1 (dynein light chain 1, LC8/PIN), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for cytoskeleton, apoptosis, and DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	DYNLL1 Rabbit mAb
Remarks/Alias	LC8; PIN; DLC1; DLC8; LC8a; DNCL1; hdlc1; DNCLC1; DYNLL1
Species	Human
Gene ID (Human)	8655
Gene ID	8655
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-89 of human DYNLL1 (P63167).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63167
Protein Weight	10kDa
Shipping	Ice bag