

Anti-Dynll2 Polyclonal Antibody For Wb And Elisa - Q96Fj2

Item Number: CM0015762



Introduction

Rabbit polyclonal antibody against human DYNLL2, validated for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for dynein and intracellular transport research.

[Learn More](#)

Parameter	Value
Product Name	DYNLL2 Rabbit pAb
Remarks/Alias	Dlc2, DNCL1B, RSPH22, DYNLL2
Species	Human
Gene ID (Human)	140735
Gene ID	140735
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96FJ2
Protein Weight	10kDa
Shipping	Ice bag