

Anti-Dyrk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92630

Item Number: CM0027064



Introduction

Anti-DYRK2 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse DYRK2, a dual-specificity kinase involved in cell cycle regulation, apoptosis, and DNA damage response. Ideal for kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	DYRK2 Rabbit pAb
Remarks/Alias	DYRK2
Species	Human
Gene ID (Human)	8445
Gene ID	8445
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 429-528 of human DYRK2 (NP_003574.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92630
Protein Weight	67kDa
Shipping	Ice bag