

Anti-Dyrk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92630

Item Number: CM0021146



Introduction

Validated DYRK2 Rabbit Polyclonal Antibody for WB, IHC-P, ELISA, cross-reacting with human and rat. Targets dual-specificity kinase involved in cell cycle, apoptosis, and DNA damage response. Research-use focused.

[Learn More](#)

Parameter	Value
Product Name	DYRK2 Rabbit pAb
Remarks/Alias	DYRK2; dual specificity tyrosine-phosphorylation-regulated kinase 2
Species	Human
GeneID	8445
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human DYRK2.
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q92630
Protein Weight	59kDa/66kDa
Shipping	Ice bag