

Anti-Dyrk4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nr20

Item Number: CM0023363

Introduction

Anti-DYRK4 rabbit polyclonal antibody validated for WB and ELISA applications. Detects dual specificity tyrosine-phosphorylation-regulated kinase 4 in human, mouse, and rat samples. Suitable for spermiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DYRK4 Rabbit pAb
Remarks/Alias	DYRK4
Species	Human
GeneID (human)	8798
GeneID	8798
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 330-480 of human DYRK4 (NP_003836.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR20
Protein Weight	60kDa
Shipping	Ice bag