


Anti-Dystrophin Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11532

Item Number: CM0016991



Introduction

Rabbit polyclonal antibody targeting human dystrophin (Gene ID: 1756) for WB, IF-P, IHC-P, ELISA. Validated in human, mouse, rat. Ideal for neuromuscular disease research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Dystrophin Rabbit pAb
Remarks/Alias	BMD; CMD3B; MRX85; DXS142; DXS164; DXS206; DXS230; DXS239; DXS268; DXS269; DXS270; DXS272; Dystrophin
Species	Human
GeneID (Human)	1756
GeneID	1756
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 346-635 of human Dystrophin (NP_003997.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11532
Protein Weight	427kDa
Shipping	Ice bag