

Anti-E-Cadherin Polyclonal Antibody For Wb, Ihc-P, Elisa - P12830

Item Number: CM0026248



Introduction

Rabbit polyclonal antibody targeting human E-Cadherin (CDH1). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse, rat. 97 kDa. Ideal for cell adhesion and cancer research.

[Learn More](#)

Parameter	Value
Product Name	E-Cadherin Rabbit pAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324; E-Cadherin
Species	Human
Gene ID	999
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 373-622 of human E-Cadherin (NP_004351.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12830
Protein Weight	97kDa
Shipping	Ice bag