

Anti-E-Cadherin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12830

Item Number: CM0014358



Introduction

Anti-E-Cadherin Rabbit pAb validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for studying cell adhesion, epithelial integrity, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	E-Cadherin Rabbit pAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324; E-Cadherin
Species	Human
GeneID (human)	999
GeneID	999
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 373-622 of human E-Cadherin (NP_004351.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12830
Protein Weight	97kDa
Shipping	Ice bag