

Anti-E-Cadherin Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P12830

Item Number: CM0002230



Introduction

High-quality rabbit monoclonal antibody against human E-Cadherin. Validated for WB, IF, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for cell adhesion, EMT, and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	E-Cadherin Rabbit mAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324; E-Cadherin
Species	Human
GeneID	999
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 701-800 of human E-Cadherin (NP_004351.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12830
Protein Weight	97kDa
Shipping	Ice bag